

The ISQ value of 1–99 is a measurement of the stability of the implant – the higher the value, the more stable the implant. If the value is 70 ISQ or above, immediate restoration is possible. In the event of values between 55 and 70, we recommend interlocked immediate restoration. If the value is below 55 ISQ, the implant must undergo submerged healing. The breident Penguin^{RFA} measures the ISQ value to an accuracy of +/- 1 ISQ unit.



Testimonial



Prof. Dr. José Eduardo Maté Sánchez de Val. PhD, MSc, DDS.

"In my opinion, the device should be used as standard, as it represents a safe and predictable solution for the user. In my scientific and clinical experience, the ISQ value is the most reliable value for determining the stability of implants. As a rule, I carry out immediate restoration of my implants with an ISQ value of 65 or above."

Literature



Worldwide, more than 700 articles have been published in peer-reviewed journals on the subject of implant stability since 1996.

We particularly recommend the publication by Professor Lars Senneryby (Resonance Frequency Analysis for Implant Stability Measurements, INTEGRATION DIAGNOSTICS UPDATE 2015;1:1-11).

www.breidentgroup.net/penguin-de

			
REF 000244GB	REF 009913GB	REF 009912GB	REF 000200GB
			
REF 009932GB	REF 009910GB	REF 000250GB	

Mistake and subject to change reserved

009936GB-20180318



breident Penguin^{RFA}



Suitable for immediate restoration?

The question of whether the fitted implant is suitable for immediate restoration often represents a challenge for the treating dentist. The torque in Ncm is used as standard when inserting the implant. The bredent Penguin^{RFA} is the ideal solution as a supplement to the ratchet.

The bredent Penguin^{RFA} is a device used to measure the stability of implants using the resonance frequency analysis (RFA). A small, magnetic measuring pin, the MultipegTM, is screwed into the implant or the abutment for this purpose and begins to vibrate without touching it. The value measured is shown as the "Implant Stability Quotient" (ISQ) and provides information regarding the appropriate treatment of the implant.


The solution: the bredent Penguin^{RFA}



The measurement method



MultiTipegTM The MultiTipegTM is available in four different sizes and makes it possible to measure primary stability in virtually all SKY[®] implant system implants.




System overview bredent Penguin^{RFA}

Instrument Kit



MultiTipegTM



Sterile cover



Penguin^{RFA} Cover 1 set (20 pieces) REF 580PENG

Ordering service

Online order form order.bredent-medical.com
 Telephone +49 7309 872-211
 Fax order +49 7309 872-24

Benefits of the bredent Penguin^{RFA}

- **Mobility**
The wireless bredent Penguin^{RFA} offers all of the freedom of mobile work. Additional patient situations can be easily checked in other treatment rooms without being restricted to the surgery unit.
- **Process security**
The bredent Penguin^{RFA} makes it possible to carry out predictable treatment, most notably in risk patients. The diagnosis device can be used to avoid complications, carry out immediate loading and reduce unnecessary costs and treatment times. As a result, process security is increased and the risk to patients is reduced.
- **Economical accessories**
The magnetic MultiTipegTM is biocompatible, can be activated without touching it and can be sterilised up to 20 times. This means that high volumes of stock can be reduced to a minimum. This saves costs and reduces the amount of administration work.

Technical data	
Electricity supply	5 V, 1 A
Charger input	100 - 240 V, 5A
Weight of the bredent Penguin ^{RFA}	100 g
The bredent Penguin ^{RFA}	is intended for long-term operation
The bredent Penguin ^{RFA}	contains NiMH batteries