

neo.lign

Full teeth

Design chart



choose natural beauty

powered by
visio.lign

Upper anterior



Upper anterior




Scale 1:1, dimensions in mm





Minor deviations in the range of tenths of a millimeter reserved

ISO 22112:2017




Lower anterior




Upper/lower posterior



Upper/lower posterior





SHADES AND DESIGNS "VARIETY WITH SYSTEM"

neo.lign A
Anterior teeth



- 13 upper moulds
- 4 lower moulds

neo.lign P
Posterior teeth
G-designs



- multifunctional occlusal design for all occlusion concepts

L-designs



- lingualised occlusal design for full dentures

Available in sizes




small large

Available in sizes



small large



Combination tables
neo.lign Full teeth

Anterior Upper	Anterior Lower	Posterior Upper / Lower G-design	Posterior Upper / Lower L-design
F41	T33	G2	L1 / L2
B43	T33 / T35	G2	L2
C43	T33 / T35	G2	L2
A44	T33 / T35	G2 / G3	L2 / L3
F44	T33 / T35	G2 / G3	L2 / L3
I45	T33 / T35	G2 / G3	L2 / L3
H46	T33 / T35	G2 / G3	L2 / L3
S46	T35 / D38	G3	L3
I47	D38	G3	L3
D48	D38	G3 / G4	L3 / L4
M48	D38	G3 / G4	L3 / L4
D49	D38 / D41	G3 / G4	L3 / L4
B51	D41	G4	L4

neo.lign Full teeth
to novo.lign Veneers
Posterior Upper / Lower

neo.lign L-design	novo.lign W-design
L1	-
L2	W3
L3	W4
L4	W5

neo.lign Design chart

Full teeth



choose natural beauty

bre^{dent}
group

000329GB-20210629
Mistake and subject to change reserved

